



BIPHASIC DEFIBRILLATOR MONITOR E-HEART

Interchangeable reusable pads



Adult and Pediatric

Flip down display Innovated system



Allows multiple viewing angles for the screen
Offers greater comfort and safety wherever the patient is.

Invasive pads



Available in different sizes:
adult and pediatric

CONTACT US

7831 NW 72ND AVENUE,
MEDLEY, FL 33166 USA

phone: +1 305 887 7552

e-mail: info@usdefib.com

THE MOST COMPACT DEFIBRILLATOR IN THE WORLD

THERMAL PRINTER



12VDC
AMBULANCE
INPUT FOR
MOBILE ICU



EASY ACCESS CONNECTORS



H 22CM X L 24CM X W 26CM



PACEMAKER



SPO2



NON
INVASIVE
BLOOD
PRESSURE
NIBP



NIBP CUFF
(AVAILABLE
IN DIFFERENT
SIZES)



Features:

- Exponential truncated Biphasic waveform
- Portable, transport handle
- 7" colored LCD – high resolution
- Flip Down Display – innovated system
- Up to 12 ECG derivations
- Ambulance input for mobile ICU
- Thermal Printer

Optional parameters:

- Pulse Oximetry – (SPO2)
- Non invasive Blood Pressure (NIBP)
- Automated defibrillation with voice and text command (AED MODE)
- External Pace Maker
- Capnography (ETCO2)
- 360 Joules

